



PRECISION AIRE TABLES & TABLETOPS

INDUSTRIAL MACHINERY & AUTOMATION





FABREEKA® ISOLATED TABLES AND TABLETOPS SIMPLE AND SMART -THE PERFECT SOLUTION FOR **QUALITY OPTIMIZATION**



High-precision equipment only provides optimal results with every possible condition met. Our systems provide maximum vibration reduction - the result of a century of experience in vibration and shock isolation.

The Precision-Aire[™] series table provides a self-leveling work surface and low frequency isolation in laboratories, cleanrooms and inspection areas where impact noise and structural vibration can decrease performance, reliability and accuracy.

A range of sizes and work surface options are available. The tables require a continuous air supply of no more than 7 bar. The work surface is isolated by pneumatic isolators. Leveling valves provide accurate, repeatable leveling of the work surface. The pneumatic isolators utilize laminar flow damping which improves settling after movement and stability for varying load conditions.


APPLICATION AREAS FOR ISOLATORS

U,I		
Active		
PAL (special)		
		_
PAT (Precision	-Aire)	
DAL (show do set	, I	
PAL (standard)	
		-
PIM		
FLII		
Rubber-metal	isolators	
:	:	
Sheet and pan	el materia	ls
I :		-

MAIN FEATURES / BENEFITS

- Fabreeka vibration-isolated tables help to provide measuring Significantly better performance due to very and analysis equipment to full performance, reliability and accuracy
- Can support loads from 100 to 1,200 kg
- Isolate work surfaces with vertical and horizontal natural frequencies natural frequencies of 1.5 Hz
- Level control and optional very quiet air compressor
- Wide range of perfectly integrated accessories

PRECISION AIRE TABLES & TABLETOPS



BENEFITS VS. SOFT-PASSIVE SYSTEM

- low natural frequencies
- Constant performance over the entire load range
- Constant performance over the entire life time
- Measuring device always aligned even in case of load changes or displacement, by automatic level control
- Adjustable pneumatic damping for optimum adaptation to the application

INDUSTRIES USING PNEUMATIC ISOLATION DEVICES

PRODUCT LINES STANDARD - OEM - CUSTOMIZED

A wide field of industries alongside medical engineering and ultra-precision measuring applications have come to use microtechnology and nanotechnology.

MACHINES

Semiconductor chip manufacturing is a very complicated process involving tens of different machines on the production floor. Most of these processes require dedicated precision equipment, chiefly high-resolution units subject to aggressive moving-stage motion.

Semiconductor-manufacturing machines are affected by vibrations from the floor or environment, or their own moving parts.

Vibrations from various sources reduce throughput and cause excessive wear.

Fabs are designed so as to keep vibration

levels within the specs of sensitive equip-

ment. Most floors in fabs do not meet the

design spec; and even if they do, equip-

ment on the fab floor increases vibration

levels over time (Colin Gordon Associates,

SEMICON-FABS

MICROSCOPES: AFM, SEM, TEM




NANO INSTRUMENTS



PRECISION MEDICAL DEVICES

Vibration levels are usually measured before the fab is equipped, but degrade by around 10 – 15 dB over the first two years with increasing amounts of production and maintenance equipment.


Production technology developing from 7 to 5 nm node pose their own challenges by increasing demands on measuring and inspection equipment.



SEMICONDUCTORS



IVF AND LIFE SCIENCE



PRECISION MEASURING EOUIPMENT



PRECISION BALANCES

Fabs that are already struggling to meet floor vibration criteria for the current technology node will have even more difficulty in the next node.

Vibration isolation platforms are needed for the equipment to function up to spec.







INTEGRATED OEM SOLUTIONS

facturers

Our standard systems are designed for easy installation with minimal technical support required.

Our technology integrates into equipment from other manu-

STANDALONE OFF-THE-

SHELE PRODUCTS



PRODUCT OVERVIEW

PAT	Precision Aire Table
ISO Tab-L	Tabletop Isolation Systems

www.fabreeka.com	

2004).



SOLUTIONS

If our off-the-shelf products are not compatible with your needs, we are happy to tailor our technology for a custom solution.

Fabreeka® offers a wide range of pneumatic vibration isolation systems. We adapt of-the-shelf systems to specific OEM specifications on technical performance or dimensional constraints towards completely bespoke solutions that meet the most exacting customer requirements.



.12

PRECISION AIRE TABLE & ISO-TAB-L

For extremely sensitive laboratories and measuring equipment the manufacturer often specifies limit values for acceptable vibration levels.


If this is not the case, we can find guide values for amplitudes of mechanical vibrations in addition in the VDI guideline VDI 2038-2.


The "Vibration Criteria" (VC lines) listed therein are frequently used internationally for the evaluation of locations for vibration-sensitive equipment.

For modern, high-precision requirements, the VC lines have been extended to include the so-called "nano lines" [Heiland et. al., 2005].

It has been shown in extensive comparisons with high-resolution SEM and TEM that the listed structural sizes can be represented if the limit values are observed.

Depending on the natural frequency and damping of a vibration isolator, disturbing frequencies are passed through in different amounts. The dimension can be calculated via the transfer function. The graph below shows the transfer function graphically for





VC-C Appropriate standard for optical microscopes to 1000x

Lithography and inspection equipment to 1000 nm detail size -----

VC-nanoD

Very difficult to reach criterion for Nanotec REM, detail size 1nm

Suitable in most instances for demanding equipment, including many electron microscopes, detail size 100-300 nm

Extremely difficult to reach

criterion for Nanotec REM,

detail size 0,2 - 0,5 nm

VC-nanoE

VC-E

<100 nm

VC-nanoEF Most strict criterion for SEM and TEM of nanotechnology in the sub-Angstrom range

Challenging criterion to

achieve for equipment of

highest precision, detail size

VC-F Criterion for extremely quiet research rooms, more than difficult to reach


VC-G Criterion for extremely quiet research rooms, more than difficult to reach

PAT Standard PAT with Pendulum air spring **ISOTabL** vertical: 2 Hz vertical: 1,8 Hz vertical: 4,5 Hz horizontal: 5 Hz horizontal: 2 Hz horizontal: 5,2 Hz

PRECISION AIRE TABLES & TABLETOPS

different table designs. With the help of this graph, it is possible to estimate the reduction of disturbing frequencies by using this product. The inverse of the transfer function is the degree of isolation.

PAT **PRECISION AIRE TABLE**



Model No.	Size (mm)	Load capacity (kg)
PAT75/75-A4*	750 x 750 x 760	370
PAT90/75-A4*	900 x 750 x 760	350
PAT120/75-A4*	1200 x 750 x 760	320
PAT150/75-A4*	1500 x 750 x 760	260

Designs: SIX with steel plate 25 mm (covered with inox laminate), G60 with granite plate 60 mm, BB with 50 mm Breadboard top

VARIOUS OPTIONS AVAILABLE

According to the customer's requirements the working surfaces are available in granite as well as stainless steel and breadboards, with thicknesses of 60 mm, 25 mm and 50 mm. The Breadboards are available in 50 mm thickness with sealed M6 holes and 25 mm pitch. Custom hole patterns and camera/laser

The extensive range of accessories for these special tables and air suspension systems include arm supports, movable shelves, casters, and front and rear protective bars.

ports for OEM applications can also be supplied.



Heavy duty version with cast iron mounting plate



Special version with top-mounted PAL air springs and 140 mm granite top

PRECISION AIRE TABLES & TABLETOPS

The protective bars are necessary for mounting movable shelves and enclosures.

If you have special requirements and need one of our special designs, require additional equipment or simply wish some technical support, please contact one of our engineers for a non-binding and cost-free consultation.



with simple top-mounted PLM air springs



Standard version with 25 mm stainless steel plate but with additional protective bar

ISO**Tab-L** VIBRATION ISOLATION SOLUTIONS TABLETOP



ISO Tab-L tabletop platforms are ideal for desktops in laboratories and Q.C. rooms, where they effectively isolate structural and floor vibration. Isolation is provided in the vertical and horizontal axes by pneumatic isolators located beneath the work surface. A protective shroud is supplied with all models which is isolated from the working surface and protects the system from accidental jarring.

	A	IFABREEKA	6	0	
G	FABREEKA	۲	0	0	
-	-	-			

Model No.	Size (mm)	Number of spring elements**	Load capacity (kg)	Weight (kg)
ISO-13*	300 x 450 x 76	3	110	20
ISO-23*	450 x 450 x 76	3	110	25
ISO-24*	450 x 450 x 76	4	145	25
ISO-33*	450 x 600 x 76	3	110	30
ISO-34*	450 x 600 x 76	4	145	30
ISO-43*	600 x 900 x 76	3	85	60
ISO-44*	600 x 900 x 76	4	130	60
ISO-53*	600 x 900 x 76	3	360	60
ISO-54*	600 x 900 x 76	4	480	60

* Designs: Granite Series G - Polished Aluminum Series A - 25 RMS Matt Brushed Stainless Steel. ** Contact Fabreeka to find out about special finishes in terms of load capacity.



PRECISION AIRE TABLES & TABLETOPS

The compact design of the ISO Tab-L enables a working height of 75 to 95 mm and will support and isolate precision instruments weighing between 22 kg and 405 kg.

Leveling of the table surface is performed via the regulator control units with air supplied from a manual pump or a dedicated source. The pressure at each isolator can be adjusted manually to accommodate changes in load distribution.



FABREEKA® THE COMPANY


Fabreeka® has been leading the international market in shock and vibration isolation since 1936. Our facilitiy at Büttelborn near Darmstadt includes our European administration as well as installation, service, quality assurance and warehousing.

Our products can be found in specialties such as measurement and laboratory technology, building services and mechanical engineering. Our in-house and field staff provide vibration measurements in the field as well as installation services, consulting and training.

This brochure describes our pneumatic isolation table and table top isolation systems product range. Don't hesitate to contact us if you have any questions or need advice on an ideal solution for your vibration issues. Our team of qualified engineers at Fabreeka® would be pleased to discuss the matter with you on the phone or meet you at your premises.

Please refer to the last page for our contact details and locations.

GLOBAL THINKING



Many companies offer products for vibration isolation and shock control. Fabreeka® demonstrably delivers sophisticated technical solutions. Understanding customer-specific requirements and matching proven products to the solution required has been one of our strengths for than hundred years.

OUR VISION

Our vision is to ensure the best results with the highest accuracy and throughput while maintaining the lowest failure rate for every item of precision equipment on the market used in research, manufacturing and measurement.

OUR MISSION

Our mission is to improve, innovate, and supply advanced technology across the world, and to bring inspiration and confidence to scientists, developers, manufacturers and consumers by ensuring a vibration-free environment for optimal results.

THE TEAM

EXPERIENCE

Our team of mechanical and civil engineers, physicists, and experts from industry share a vast wealth of experience in automation, electronics and real-time processing.


CUSTOMIZABLE SOLUTIONS

Tailored to customer and industry needs

PATENTED TOP-PERFORMING SYSTEMS

DESIGN-FOCUSED R&D In-house manufacturing

Innovative, user-friendly



PRECISION AIRE TABLES & TABLETOPS

YOUR MOTION. OUR SOLUTION. **STABILUS**

THE GLOBAL MOTION CONTROL SOLUTION **PROVIDER TO THE INDUSTRY**

- Specialized industrial brands with unique solutions and expertise for each market.
- Wide product range from simple to complex requirements across all industries.
- Broad global presence that closely connects us with all markets.

FABREEKA[®] STABILUS ACE HAHN Globally leading series gas Leading supplier for individual Expert in industrial damping Experienced provider of spring manufacturer meeting vibration isolation and shock and deceleration with excellent customized and special gas highest automotive quality control solutions. customer support satisfying springs with **very fast** time to solutions. standards for industrial most demanding application market. products in all plants. requirements.

Stabilus group of industrial brands







Fabreeka GmbH Deutschland

Hessenring 13 64572 Büttelborn Germany T +49 6152 - 9597-0 F +49 6152 - 9597-40

info@fabreeka.de www.fabreeka.de

QR code takes you to our International Locations



in @fabreeka_international_inc f @facebook.com/FabreekaIntl